

## Product datasheet for **QP100016**

### SensiFast SYBR Master Mix Lo-ROX (2x)

#### Product data:

<b>Product Type:</b>	qPCR Reagents
<b>Components:</b>	One vial containing 1.0 mL SensiFAST SYBR Master Mix - Lo-ROX Kit (sufficient for 100 rxns) One Datasheet
<b>Quality Control:</b>	The SensiFAST SYBR Lo-ROX Kit has been developed for fast, highly accurate real-time PCR and has been validated on all commonly-used real-time instruments that require a low concentration of the passive reference dye ROX. A combination of the latest advances in buffer chemistry and PCR enhancers ensures that the SensiFAST SYBR Lo-ROX Kit produces reliable assay results under fast thermal cycling conditions. An antibody-mediated hot-start DNA polymerase system promotes highly-specific amplification, in turn improving assay sensitivity and dynamic range.
<b>Storage:</b>	Upon receipt, please store the Master Mix at -20°C.
<b>Stability:</b>	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)