

Product datasheet for **QP100019**

SensiFAST SYBR Master Mix Fluorescein (2x)

Product data:

Product Type: qPCR Reagents

Components: One vial containing 1.0 mL SensiFAST SYBR Master Mix - Fluorescein kit (sufficient for 100 rxns) One Datasheet

Quality Control: The SensiFAST SYBR & Fluorescein Kit has been developed for fast, highly accurate real-time PCR and has been validated on real-time PCR platforms that require the passive reference dye fluorescein for collection of dynamic well factors. A combination of the latest advances in buffer chemistry and PCR enhancers ensures that the SensiFAST SYBR & Fluorescein Kit produces reliable assay results under fast thermal cycling conditions. An antibody-mediated hot-start DNA polymerase system promotes highly-specific amplification, in turn improving assay sensitivity and dynamic range.

Storage: Upon receipt, please store the Master Mix at -20°C.

Stability: The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)