

## Product datasheet for **QP100024**

### SensiFAST SYBR Master Mix No-ROX (2x)

#### Product data:

<b>Product Type:</b>	qPCR Reagents
<b>Components:</b>	Ten vials of QP100022S each containing 1.0 mL SensiFAST SYBR Master Mix - No-ROX Kit (sufficient for 1000 rxns) One Datasheet
<b>Quality Control:</b>	The SensiFAST SYBR No-ROX Kit has been developed for fast, highly accurate real-time PCR and has been validated on all commonly-used real-time instruments that do not require the passive reference dye ROX. A combination of the latest advances in buffer chemistry and PCR enhancers ensures that the SensiFAST SYBR No-ROX Kit produces reliable assay results under fast thermal cycling conditions. An antibody-mediated hot-start DNA polymerase system promotes highly-specific amplification, in turn improving assay sensitivity and dynamic range.
<b>Storage:</b>	Upon receipt, please store the Master Mix at -20°C.
<b>Stability:</b>	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)