

Product datasheet for R1020

RFX5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, EMSA, IP
Recommended Dilution:	Suitable for ELISA and supershift assays. This product was also tested in a gel supershift assay and found to be reactive against RFX5 complexes using 0.5 to 1.0 ml per assay.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	RFX5 peptide corresponding to amino acids 320 to 494 of the human protein conjugated to Keyhole Limpet Hemocyanin (KLH).
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with 0.01% sodium azide as preservative. State: Serum State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Prepared from monospecific antiserum by delipidation and defibrination.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid cycles of freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	regulatory factor X5
Database Link:	Entrez Gene 5993 Human P48382
Synonyms:	Regulatory factor X 5



[View online »](#)