

Product datasheet for R1104

Penicillinase Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Suitable for Immunoblotting (Western or Dot blot), ELISA and most Immunological methods requiring high titer and specificity. <u>Recommended Dilutions:</u> This product has been assayed against 1.0 ug of Penicillinase [Enterobacter cloacae] in a standard sandwich ELISA using Peroxidase conjugated Affinity Purified anti- Rabbit IgG [H&L] (Goat) (R1364HRP) and (ABTS (2,2'-azino-bis-[3-ethylbenthiiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature. A working dilution of 1:6,000 to 1:26,000 of the reconstitution concentration is suggested for this product.
Reactivity:	Enterobacter
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Penicillinase [Enterobacter cloacae].
Specificity:	Assay by immunoelectrophoresis resulted in a single precipitin arc against purified and partially purified Penicillinase [Enterobacter cloacae]. Cross reactivity against Penicillinase from other tissues and species may occur but have not been specifically determined.
Formulation:	0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4 with 0.01% (w/v) sodium azide as preservative. State: Serum State: Lyophilized purified Ig fraction.
Reconstitution Method:	Restore with 2.0 ml of deionized water (or equivalent).
Concentration:	lot specific
Purification:	Prepared from monospecific antiserum by a delipidation and defibrination.
Conjugation:	Unconjugated



[View online »](#)

Storage: Store vial at 2-8°C prior to restoration. For extended storage add glycerol to 50% and then aliquot contents and freeze at -20°C or below. Centrifuge product if not completely clear after standing at room temperature.
This antibody is stable for one month at 2-8°C as an undiluted liquid.
Dilute only prior to immediate use.
Avoid repeated freezing and thawing.

Stability: Shelf life: One year from despatch.