

Product datasheet for **R1114HRPS**

KT13 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Western blot: 1/1,000-1/5,000. ELISA: 1/10,000-1/40,000. Antibody has been assayed against 1.0 µg of Trypsin inhibitor [lima bean] in a standard capture ELISA using ABTS as a substrate for 30 minutes at room temperature. A working dilution of 1/20,000 to 1/100,000 of the reconstitution concentration is suggested.
Reactivity:	Lima Bean
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Trypsin inhibitor from lima bean
Specificity:	This antibody detects lima bean Trypsin inhibitor. Cross reactivity against Trypsin inhibitor from other sources may occur but have not been specifically determined. Immunoelectrophoresis give a single precipitin arc against anti-peroxidase and anti-rabbit serum as well as purified and partially purified Trypsin inhibitor [lima bean].
Formulation:	0.02 M Potassium phosphate, 0.15 M Sodium chloride, pH 7.2 Label: HRP State: Lyophilized IgG fraction Stabilizer: 10 mg/ml BSA (immunoglobulin and protease free) Preservative: 0.01% (w/v) Gentamicin sulfate (Do NOT add Sodium azide!) Label: Horseradish peroxidase
Reconstitution Method:	Restore with 0.1 ml of deionized water or equivalent.
Concentration:	lot specific
Purification:	Delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer
Conjugation:	HRP
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.



[View online »](#)

Stability: Shelf life: one year from despatch.

Database Link: [P01070](#)