

Product datasheet for **R1150BS**

ycjM Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IP, WB
Recommended Dilution:	Western blot: 1/500-1/5,000. Immunoprecipitation: 1/100. ELISA: 1/5,000-1/20,000. This product has been assayed against 1.0 µg of Sucrose Phosphorylase [E.coli] in a standard capture ELISA using peroxidase conjugated streptavidin and ABTS as a substrate for 30 minutes at room temperature. A working dilution of 1/25,500 to 1/115,000 of the reconstitution concentration is suggested for this product.
Reactivity:	Escherichia coli
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Sucrose phosphorylase from E.coli
Specificity:	This antibody detects E.coli Sucrose phosphorylase. Cross reactivity against Sucrose phosphorylase from other sources is unknown. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-goat serum as well as purified and partially purified Sucrose phosphorylase [E.coli].
Formulation:	0.02 M Potassium phosphate, 0.15 M Sodium chloride, pH 7.2 Label: Biotin State: Purified State: Lyophilized purified Ig fraction Stabilizer: 10 mg/m BSA (immunoglobulin and protease free) Preservative: 0.01% (w/v) Sodium azide
Reconstitution Method:	Restore with 0.1 ml of deionized water (or equivalent).
Concentration:	lot specific
Purification:	Delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer
Conjugation:	Biotin



[View online »](#)

Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P76041
Synonyms:	Sucrose glucosyltransferase (E. coli)