

Product datasheet for **R1250P**

Monkey IgA + IgG + IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Monkey IgA + IgG + IgM Goat Polyclonal Antibody
Recommended Dilution:	Immunoprecipitation. Immunodiffusion. ELISA: 1/100000. Western Blot: 1/2000 - 1/10000.
Host:	Goat
Immunogen:	Monkey IgG, IgA and IgM whole molecules
Formulation:	0.125 M Sodium Borate, 0.075 M Sodium Chloride, 0.005 M EDTA; pH 8.0 containing 0.01% Sodium Azide State: Aff - Purified State: Liquid
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. For extended storage add glycerol to 50%. Avoid repeated freezing and thawing.



[View online »](#)