

## Product datasheet for **R1255AP**

### Rabbit IgG (H+L chain) Sheep Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG (H+L chain) Sheep Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoblotting (Western or Dot blot), ELISA and Immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency. <u>Recommended Dilutions:</u> ELISA: 1/15,000. Western Blot: 1/500-1/2,000. Immunohistochemistry: 1/200-1/1,000.
Reactivity:	Rabbit
Host:	Sheep
Immunogen:	Rabbit IgG whole molecule.
Formulation:	0.05M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol, pH 8.0 with 0.01% Sodium Azide as preservative and 10 mg/ml BSA (IgG and Protease free) as stabilizer. Label: AP State: Liquid (sterile filtered) purified Ig fraction. Label: Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
Concentration:	lot specific
Purification:	Prepared from monospecific antiserum by Immunoaffinity Chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity. Dilute only prior to immediate use.



[View online »](#)