

Product datasheet for **R1302P**

Chicken IgG / Chicken IgY (F(ab)₂ specific) Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (F(ab) ₂ specific) Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, Conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity.
Reactivity:	Chicken
Host:	Rabbit
Immunogen:	Chicken IgG F(ab') ₂ fragment.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with 0.01% (w/v) sodium azide as preservative. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)