

Product datasheet for **R1312P**

Canine IgG (F(ab)2 specific) Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgG (F(ab)2 specific) Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, Conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. <u>Recommended Dilutions:</u> ELISA: 1/20,000-1/100,000 Western blot: 1/2,000-1/10,000. Immunohistochemistry: 1/1,000-1/5,000.
Reactivity:	Dog
Host:	Rabbit
Immunogen:	Dog IgG F(ab') ₂ fragment.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 without preservatives. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)