

Product datasheet for R1312T

Canine IgG (F(ab)2 specific) Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine lgG (F(ab)2 specific) Rabbit Polyclonal Antibody
Recommended Dilution:	FLISA: 1/10000 - 1/50000. Immunofluorescence: 1/1000 - 1/5000. Flow Cytometry: 1/500 - 1/2500.
Reactivity:	Dog
Host:	Rabbit
Immunogen:	Dog IgG F(ab')2 fragment
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 Label: TRITC State: Lyophilized Canine IgG F(ab')2 Label: Rhodamine Absorption emission: 550 nm / 570 nm
Reconstitution Method:	Restore with 1.0 mL of deionized water (or equivalent).
Concentration:	lot specific
Conjugation:	TRITC
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US