

Product datasheet for **R1325P**

Guinea pig IgG (H+L chain), adsorbed Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Guinea pig IgG (H+L chain), adsorbed Goat Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, Conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. A dilution of 1:50,000 is suggested for ELISA.
Reactivity:	Guinea Pig
Host:	Goat
Immunogen:	Guinea Pig IgG whole molecule.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with 0.01% sodium azide as preservative. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)