

## Product datasheet for **R1331HRP**

### Equine IgG (F(ab)2 specific) Rabbit Polyclonal Antibody

#### Product data:

<b>Product Type:</b>	Secondary Antibodies
<b>Product Name:</b>	Equine IgG (F(ab)2 specific) Rabbit Polyclonal Antibody
<b>Recommended Dilution:</b>	Suitable for Immunoblotting (Western or Dot blot: use at 1/1,000-1/5,000), ELISA, Immunoperoxidase electron microscopy and Immunohistochemistry (use at 1/500-1/2,500) as well as other peroxidase-antibody based enzymatic assays requiring lot-to-lot consistency. This product has been assayed against 1.0 µg of Horse IgG in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature. A working dilution of 1:10,000 to 1:50,000 of the reconstitution concentration is suggested for this product.
<b>Reactivity:</b>	Horse
<b>Host:</b>	Rabbit
<b>Immunogen:</b>	Horse IgG F(ab')2 fragment.
<b>Formulation:</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with 10 mg/ml Bovine Serum Albumin (BSA, IgG and Protease free) as stabilizer and 0.01% (w/v) Gentamicin Sulfate as preservative Label: HRP State: Lyophilized purified Ig fraction Label: Horseradish Peroxidase
<b>Reconstitution Method:</b>	Restore with 1.0 ml of deionized water (or equivalent).
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Immunoaffinity Chromatography
<b>Conjugation:</b>	HRP
<b>Storage:</b>	Store vial at 2-8°C prior to restoration. For extended storage add glycerol to 50% and then aliquot contents and freeze at -20°C or below. Centrifuge product if not completely clear after standing at room temperature. This antibody is stable for one month at 2-8°C as an undiluted liquid. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)