

Product datasheet for **R1337P**

Human IgG (Fc specific) Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Rabbit Polyclonal Antibody
Applications:	ELISA, IHC, IP, WB
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, Conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. <u>Recommended dilutions:</u> ELISA: 1:10,000-1:50,000. Western Blot: 1:1,000-1:5,000. Immunohistochemistry: 1:500-1:3,000.
Reactivity:	Human
Host:	Rabbit
Immunogen:	Human IgG F(c) fragment
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide as preservative State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)