

Product datasheet for **R1339P**

Human IgG (F(ab)2 specific), adsorbed Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (F(ab)2 specific), adsorbed Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, Conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. <u>Recommended Dilutions:</u> ELISA: 1/350,000. Western Blot: 1/3,000-1/15,000. Immunohistochemistry: 1/1,000-1/5,000.
Reactivity:	Human
Host:	Rabbit
Immunogen:	Human IgG F(ab')2 fragment.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide as preservative. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)