

Product datasheet for R1342AP

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

Human IgA Rabbit Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Human IgA Rabbit Polyclonal Antibody

Recommended Dilution: Suitable for Immunoblotting (Western or Dot blot), ELISA and Immunohistochemistry as well

as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency. Recommended Dilutions: This product has been assayed against 1.0 ug of Human IgA in a standard capture ELISA using p-nitrophenyl phosphate (pNPP) as a substrate for 30 minutes

at room temperature. A working dilution of 1:1,000 to 1:3,000 of the reconstitution

concentration is suggested for this product.

Reactivity: Human Host: Rabbit

Immunogen: Human IgA heavy chain.

Formulation: 0.05M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc

Chloride, 50% (v/v) Glycerol, pH 8.0 with 10 mg/ml BSA (IgG and Protease free) as stabilizer

and 0.01% (w/v) Sodium Azide as preservative.

Label: AP

State: Liquid (sterile filtered) purified IgG fraction.

Label: Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)

Concentration: lot specific

Purification: Immunoaffinity Chromatography.

Conjugation: AP

Storage: Store the antibody (undiluted) at 2-8°C.

DO NOT FREEZE!

Freezing alkaline phosphatase conjugates will result in a substantial

loss of enzymatic activity.

Dilute only prior to immediate use.

