

Product datasheet for **R1344P**

Monkey IgG (gamma chain specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Monkey IgG (gamma chain specific) Goat Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. ELISA: 1:20,000 - 1:100,000 dilution. Western Blot: 1:2000 - 1:10,000 Immunohistochemistry: 1:1000-1:5000
Reactivity:	Monkey
Host:	Goat
Immunogen:	Monkey IgG γ chain.
Formulation:	0.02M potassium phosphate, 0.15M sodium chloride, pH 7.2, 0.01% sodium azide. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)