

Product datasheet for **R1345AP**

Monkey IgM (Mu chain specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Monkey IgM (Mu chain specific) Goat Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoblotting (1/500-1/3,000; Western or Dot blot), ELISA (1/2,000-1/10,000) and Immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency.
Reactivity:	Monkey
Host:	Goat
Immunogen:	Monkey IgM μ chain.
Formulation:	0.05M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, pH 8.0, with 10 mg/ml Bovine Serum Albumin (BSA, IgG and Protease free) as stabilizer and 0.01% (w/v) sodium azide as preservative. Label: AP State: Liquid (sterile filtered) purified Ig fraction. Label: Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	AP
Storage:	Store the antibody (undiluted) at 2-8°C. DO NOT FREEZE! Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity. Dilute only prior to immediate use.



[View online »](#)