

## Product datasheet for R1354AP

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Mouse IgG (F(ab)2 specific), adsorbed Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Secondary Antibodies

**Product Name:** Mouse IgG (F(ab)2 specific), adsorbed Rabbit Polyclonal Antibody

Recommended Dilution: Suitable for Immunoblotting (Western or Dot blot), ELISA and Immunohistochemistry as well

as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency. Recommended Dilutions: This product has been assayed against 1.0 µg of Mouse IgG in a standard capture ELISA using pNPP p-nitrophenyl phosphate as a substrate for 30 minutes at room temperature. A working dilution of 1:500 to 1:3,500 of the reconstitution concentration

is suggested for this product.

Reactivity: Mouse Host: Rabbit

**Immunogen:** Mouse IgG F(ab')2 fragment.

**Formulation:** 0.05M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc

Chloride, 50% (v/v) Glycerol; pH 8.0 with 10 mg/ml Bovine Serum Albumin (BSA) IgG and

Protease free as stabilizer and 0.01% (w/v) sodium azide as preservative.

Label: AP

State: Liquid (sterile filtered) purified Ig fraction.

Label: Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)

**Concentration:** lot specific

**Purification:** Immunoaffinity chromatography.

Conjugation: AP

Storage: Store the antibody (undiluted) at 2-8°C.

DO NOT FREEZE!

Freezing alkaline phosphatase conjugates will result in a substantial

loss of enzymatic activity.

Dilute only prior to immediate use.

