

## Product datasheet for **R1360P**

### Mouse IgG2a (gamma 2a chain specific) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (gamma 2a chain specific) Rabbit Polyclonal Antibody
Recommended Dilution:	This antibody is suitable for immunoprecipitation, Immunodiffusion, ELISA (1/40,000-1/200,000), Western blot (1/4,000-1/20,000), Immunohistochemistry (1/2,000-1/10,000), Immunomicroscopy as well as other antibody based assay using streptavidin or avidin conjugates requiring lot-to-lot consistency.
Reactivity:	Mouse
Host:	Rabbit
Immunogen:	Mouse IgG2a $\gamma$ 2a chain
Formulation:	0.2 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, containing 10 mg/ml Bovine Serum Albumin (BSA, IgG and Protease free) as stabilizer, 0.01% (w/v) Sodium Azide as preservative State: Aff - Purified State: Liquid (sterile filtered) purified IgG fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)