

Product datasheet for **R1361P**

Mouse IgG2b (gamma 2b chain specific) Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2b (gamma 2b chain specific) Rabbit Polyclonal Antibody
Recommended Dilution:	Immunoprecipitation. ELISA: 1/20000 - 1/100000. Western Blot: 1/2000 - 1/10000.
Reactivity:	Mouse
Host:	Rabbit
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide State: Aff - Purified State: Liquid
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)