

## Product datasheet for **R1371P**

### Rat IgG (H+L chain) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG (H+L chain) Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. <i>Recommended Dilutions:</i> ELISA: 1/20,000-1/100,000. Western blot: 1/2,000-1/10,000. Immunohistochemistry: 1/1,000-1/5,000.
Reactivity:	Rat
Host:	Rabbit
Immunogen:	Rat IgG whole molecule.
Formulation:	0.02M Potassium Phosphate, 0.15M Sodium Chloride, pH 7.2 State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction. Stabilizer: None Preservative: 0.01% (w/v) Sodium Azide
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid cycles of freezing and thawing.



[View online »](#)