

Product datasheet for **R1376P**

Rat IgA (alpha chain) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgA (alpha chain) Goat Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity. <u>Recommended Dilutions:</u> ELISA: 1:10,000-1:50,000. Western Blot: 1:1,000-1:5,000. Immunohistochemistry: 1:500-1:3,000.
Reactivity:	Rat
Host:	Goat
Immunogen:	Rat IgA α chain
Formulation:	0.1 M Sodium Borate, 0.005 M EDTA, 0.25 M Sodium Chloride, pH 8.2 containing 0.01% (w/v) Sodium Azide as preservative. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid cycles of freezing and thawing.



[View online »](#)