

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for R1377F

Rat IgA + IgG + IgM (H+L chain) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat lgA + lgG + lgM (H+L chain) Goat Polyclonal Antibody
Recommended Dilution:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency. <i>Recommended Dilutions</i> : FLISA: 1/10,000-1/50,000. Flow Cytometry: 1/2,000-1/10,000. Immunofluorescence microscopy: 1/1,000-1/5,000.
Reactivity:	Rat
Host:	Goat
Immunogen:	Rat IgG, IgA and IgM whole molecules.
Formulation:	 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 Label: FITC State: Lyophilized purified Ig fraction. Stabilizer: 10 mg/ml BSA (IgG and Protease free) Preservative: 0.01% (w/v) Sodium Azide Label: Fluorescein isothiocyanate (Molecular Weight 390 daltons) Absorption emission: 495 nm / 528 nm Molar radio: 3.5 moles FITC per mole of Goat IgG.
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent).
Concentration:	2.0 mg/ml (by UV absorbance at 280 nm)
Purification:	Immunoaffinity chromatography.
Conjugation:	FITC
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US