

Product datasheet for **R1383P**

Sheep IgG F(ab)2 specific (adsorbed) Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Sheep IgG F(ab)2 specific (adsorbed) Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity.
Reactivity:	Sheep
Host:	Rabbit
Immunogen:	Sheep IgG F(ab')2 fragment.
Formulation:	in 0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.6 with 0.01% (w/v) Sodium Azide as preservative and no stabilizers. State: Aff - Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid cycles of freezing and thawing.



[View online »](#)