

Product datasheet for **R1385P**

Porcine IgG F(c) specific Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG F(c) specific Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for Immunoprecipitation, Immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity.
Host:	Rabbit
Immunogen:	Swine IgG F(c) fragment.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide as preservative State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Dilute only prior to immediate use.



[View online »](#)