

Product datasheet for R1388T

Human IgG (H+L chain), Fab Fragment Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (H+L chain), Fab Fragment Goat Polyclonal Antibody
Recommended Dilution:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Reactivity:	Human
Host:	Goat
Immunogen:	Human IgG whole molecule
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, containing 10 mg/ml BSA, IgG and Protease free, and 0.01% (w/v) Sodium Azide Label: TRITC State: Lyophilized Fab fragments Label: Tetramethylrhodamine isothiocyanate (Molecular Weight 444 daltons) Absorption emission: 550 nm / 570 nm Molar ratio: 1.2 moles TRITC per mole of goat IgG Fab
Reconstitution Method:	Restore with 1.0 ml of deionized water or equivalent.
Purification:	Immunoaffinity chromatography
Conjugation:	TRITC
Storage:	Store vial at 2-8°C prior to restoration. Restore with deionized water (or equivalent); centrifuge product if not completely clear after standing at room temperature. This product is stable for one month at 2-8°C as an undiluted liquid. For extended storage aliquot contents and freeze at -20°C or below. Avoid repeated freezing and thawing. Dilute only prior to immediate use.



[View online »](#)