

Product datasheet for **R1413P**

Rat IgG (H+L chain), F(ab)2 Fragment Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG (H+L chain), F(ab)2 Fragment Rabbit Polyclonal Antibody
Recommended Dilution:	Suitable for highly specific immunological methods requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. <u>Recommend Dilutions:</u> ELISA: 1/20,000-1/100,000 Western Blot: 1/2,000-1/10,000 Immunohistochemistry: 1/1,000-1/5,000.
Reactivity:	Rat
Host:	Rabbit
Immunogen:	Rat IgG whole molecule.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide as preservative. State: Aff - Purified State: Liquid (sterile filtered) F(ab')2 fragments.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (undiluted) at 2-8°C for one month or (in aliquots) at -20°C for longer. Dilute only prior to immediate use. Avoid repeated freezing and thawing.



[View online »](#)