

## Product datasheet for **R1456TR**

### Sheep IgG (H+L chain) Rabbit Polyclonal Antibody

#### Product data:

<b>Product Type:</b>	Secondary Antibodies
<b>Product Name:</b>	Sheep IgG (H+L chain) Rabbit Polyclonal Antibody
<b>Recommended Dilution:</b>	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency.
<b>Reactivity:</b>	Sheep
<b>Host:</b>	Rabbit
<b>Immunogen:</b>	Sheep IgG whole molecule
<b>Formulation:</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, with 10 mg/ml BSA (IgG and Protease free) and 0.01% (w/v) Sodium Azide Label: Texas Red State: Lyophilized IgG fraction Label: (TM) Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption emission: 596 nm / 620 nm Molar ratio: 2.3 moles Texas Red (TM) per mole of Rabbit IgG
<b>Reconstitution Method:</b>	Restore with 2.0 ml of deionized water or equivalent.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Immunoaffinity chromatography
<b>Conjugation:</b>	Texas Red
<b>Storage:</b>	Store vial at 4°C prior to restoration. Restore with deionized water (or equivalent); centrifuge product if not completely clear after standing at room temperature. This product is stable for four weeks at 4°C as an undiluted liquid. For extended storage reconstitute product with 50% glycerol instead of water and then aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Dilute only prior to immediate use.



[View online »](#)