

Product datasheet for R1481P

RAD9A pSer1260 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	<p>This phospho specific polyclonal antibody was tested by Immunoblotting and ELISA. Data from both Immunoblotting and ELISA indicate the antibody is reactive with the phosphorylated form of the immunizing peptide and minimally reactive with the non-phosphorylated form of the immunizing peptide. Immunoblotting detects yeast Rad9 protein. No reactivity is expected against the human or mouse analogs of RAD9. Reactivity against RAD9 from other sources is unknown. Cross reactivity may occur with auto-phosphorylated Rad53 kinase. Although not tested, this antibody is likely functional by IHC and IP.</p> <p><u>Recommended Dilution(s):</u> This product has been assayed against 0.1 µg of phosphorylated peptide (pS1260) in a standard capture ELISA using TMB (3,3',5,5'-Tetramethylbenzidine) as a substrate for 30 minutes at room temperature. A working dilution of 1:5,000 is suggested for this product. Minimal reactivity was detected against the nonphosphorylated form (S1260) of the immunizing peptide. This antibody appears to be specific for the active form (phosphorylated) of the protein. Dilute the antibody 1:100 to 1:500 for immunoblotting.</p>
Reactivity:	Yeast
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to phosphorylated form of aa 1249-1263 of 1309 of yeast Rad9 protein conjugated to KLH.
Specificity:	This affinity purified antibody is directed against the phosphorylated form of yeast Rad9 at the pS1260 residue. The product was affinity purified from monospecific antiserum by immunoaffinity purification. Antiserum was first purified against the phosphorylated form of the immunizing peptide. The resultant affinity purified antibody was then cross-adsorbed against the non-phosphorylated form of the immunizing peptide. This phospho specific polyclonal antibody reacts with phosphorylated pS1260 of yeast Rad9. Reactivity with non-phosphorylated yeast Rad9 is minimal by ELISA and immunoblotting. No reactivity is expected against the human or mouse analogs of RAD9. Reactivity against RAD9 from other sources is unknown. Cross reactivity may occur with auto-phosphorylated Rad53 kinase.



[View online »](#)

Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with 0.01% (w/v) Sodium Azide as preservative. State: Aff - Purified State: Liquid (sterile filtered) purified Ig fraction.
Concentration:	lot specific
Purification:	Immunoaffinity chromatography.
Conjugation:	Unconjugated
Storage:	Store vial at -20°C. For extended storage aliquot contents and freeze at -20°C or below. Dilute only prior to immediate use. Avoid cycles of freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	RAD9 checkpoint clamp component A
Database Link:	Q99638
Background:	<p>Rad9 is required for the MEC1/TEL1-dependent activation of <i>Saccharomyces cerevisiae</i> DNA damage checkpoint pathways mediated by Rad53 and Chk1. DNA damage induces Rad9 phosphorylation, and Rad53 specifically associates with phosphorylated Rad9. Cells have evolved multiple strategies for tolerating genomic damage. The most important of these are numerous repair systems that remove or bypass potentially mutagenic DNA lesions. Another cellular strategy is to delay cell-cycle transitions at multiple points. The genetic control of these delays, termed 'checkpoints', was first established in budding yeast where it was shown that the RAD9 gene functions in G2/M arrest after irradiation with X-rays. Subsequently, it has become clear that Rad9 also functions at the G1/S, intra-S and mid-anaphase checkpoints.</p> <p>Defects in checkpoint regulation can lead to genome instability and, in higher eukaryotes, neoplastic transformation. Rad9 also controls the transcriptional induction of a DNA damage regulon (DDR). Rad9 may also have a pro-apoptotic function. This is suggested in that Rad9 from <i>Schizosaccharomyces pombe</i> (SpRad9) contains a group of amino acids with similarity to the Bcl-2 homology 3 death domain, which is required for SpRad9 interaction with human Bcl-2 and apoptosis induction in human cells.</p> <p>Overexpression of Bcl-2 in <i>S. pombe</i> inhibits cell growth independently of rad9, but enhances resistance of rad9-null cells to methyl methanesulfonate, ultraviolet and ionizing radiation. Rad9 conveys the checkpoint signal by activating Rad53p and Chk1p; is hyperphosphorylated by Mec1p and Tel1p; and is a potential Cdc28p substrate. Mature yeast Rad9 is reported to have an apparent molecular weight of ~148 kDa. The human homolog is reported at 48.5 kDa.</p>
Synonyms:	RAD-9A, Cell cycle checkpoint control protein RAD9A, EC=3.1.11.2, DNA repair exonuclease rad9 homolog A

Note:
Protein Sequence: *Saccharomyces cerevisiae*

1 msgqlvqwks spdrvtqsai kealhsplad gdmnemnvpv dplenkvnst niiegspkan
 61 pnpvkfmnts eifqkslgll desprhddel nievgdndrp nanilhnert pdldrianff
 121 ksnrtpgken lltkyqssdl edtplmrkk mtfqtptdpl eqktfkklls dtgfcyygeq
 181 ndgeenasle vteadatfvq maersadnyd calegivtpk rykdelsksg gmqdervqkt
 241 qimisaespn sissydknki tgngrtrrnv nkvfnnnedn igaieeknpv kkkxenysd
 301 dlrrennqii qsneseeine leknlvsvgr endvnnldid insavsgtps rnaaeemys
 361 sesvnnreps kkwifryskd ktensnrst qivnnpqtqe mpldsisidt qplsksfnte
 421 tneletqii vsslsqgisa qkgpvfhtsg qteeiktqii nspeqalna tfetpvtlrs
 481 infepilevp etsspskntm skpsnsspip kekdtfnihe revetnnvfs ndiqnssnaa
 541 trddiiiags sdfneqkeit driylqlsgk qisdsgsdet ermsspeldt kkestimsev
 601 eltqelpeve eqqdlqtspk klvveetlm eikkskgnsl qlhddnkecn sdkqdgtesl
 661 dvaliehesk gqsselqknl mqlfposesqe iiqnrtrtkr rqkdtieige eenrstkts
 721 ptkhlkrnsd ldaasikrep scsitiqtge tsggkdskeq syvfpegirt adnsflskdd
 781 iifgnavwcq ytwnykfypg illevdtnqd gcwiyfetgr sltkdediyy ldirigdavt
 841 fdgneyvvvg lecrshdni ircirgydtv hlkkknasgl lgkrtlikal ssislsew
 901 akrakiiled neknkgdayr ylrhpirgrk smtnvlspkk htdekdint htevynneie
 961 sssekkeivk kdsrdalaeh agapsllfss geirtgnvfd kcifvltslf enreelrqi
 1021 esqggvties gfstlfnfth plakslvng ntdnirelal klawkphslf adcrfacit
 1081 krlhrlskyl etlalgwptl hwkfisacie kkrivphliy qyllpsgesf rlsldspsgk
 1141 giiksnnifs fytqflrgsn lrdqicgvkk mldyivivw grseldsvk fafaclsagr
 1201 mltidlpid vddtepllna ldslvprigs elsnrklkfl iyanennngks qmkllrlrs
 1261 qislkfkknf yifhteskew liqtiinedt gfhdditdnd iyntisevr


Product images:


Figure 1. Immunoblotting. Affinity purified phosphospecific antibody to yeast Rad9 at pS1260 was used at a 1:200 dilution incubated overnight at 4C to detect Rad9 by Western blot. Lanes were loaded with 50 ng each of recombinant GST fusion protein containi

Schematic summary of the DNA replication and DNA damage checkpoints in *S. cerevisiae*.




Figure 2. Summary Table. Checkpoints are mechanisms that impose delays in the cell cycle in response to DNA damage or defects in DNA replication, to ensure that mitotic transmission is error-free. Failure to delay the cell cycle in the presence of damage c