

Product datasheet for **RA020AG**

Chicken IgG / Chicken IgY Chicken Protein

Product data:

Product Type:	Native Proteins
Description:	Chicken IgG / Chicken IgY chicken protein, 10 mg
Species:	Chicken
Protein Source:	Serum
Purity:	>95%
Conjugation:	Agarose
Buffer:	State: Liquid Suspension of agarose beads Buffer System: 0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4 containing 0.01% (w/v) Sodium Azide Presentation Label: Agarose
Preparation:	Liquid Suspension of agarose beads
Applications:	Immunoprecipitation.
Protein Description:	This product is an purified Chicken IgG isolated from normal Chicken Serum coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of anti-Chicken IgG.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)