

Product datasheet for **RA020AP**

Chicken IgG / Chicken IgY Chicken Protein

Product data:

Product Type:	Native Proteins
Description:	Chicken IgG / Chicken IgY chicken protein, 1 mg
Species:	Chicken
Protein Source:	Serum
Purity:	>95% prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by conjugation and extensive dialysis against the buffer stated above
Conjugation:	AP
Buffer:	State: Liquid - sterile filtered serum Buffer System: 0.05M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol; pH 8.0 <u>Label:</u> Alkaline Phosphatase (Calf Intestine) Presentation Label: AP
Preparation:	Liquid - sterile filtered serum
Protein Description:	Alkaline Phosphatase conjugated Chicken IgG whole molecule.
Storage:	Store the product undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)