

## Product datasheet for **RA022AP**

### Bovine IgG Bovine Protein

#### Product data:

Product Type:	Native Proteins
Description:	Bovine IgG bovine protein, 1 mg
Species:	Bovine
Protein Source:	Serum
Purity:	>95% delipidation, salt fractionation, ion exchange chromatography followed by conjugation and extensive dialysis against the buffer stated above
Conjugation:	AP
Buffer:	State: Liquid Buffer System: 0.05M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol; pH 8.0 <u>Label:</u> Alkaline Phosphatase (Calf Intestine) Presentation Label: AP
Preparation:	Liquid
Protein Description:	Bovine IgG
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)