

Product datasheet for **RA023AG**

Feline IgG Feline Protein

Product data:

| | |
|----------------------|--|
| Product Type: | Native Proteins |
| Description: | Feline IgG feline protein, 10 mg |
| Species: | Feline |
| Protein Source: | Serum |
| Purity: | >95% |
| Conjugation: | Agarose |
| Buffer: | State: Liquid Buffer System: 0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4 containing 0.01% (w/v) Sodium Azide Presentation Label: Agarose |
| Preparation: | Liquid |
| Applications: | Suitable for antibody purification or to remove unwanted cross-reactivity. |
| Protein Description: | Feline IgG |
| Storage: | Store the antibody undiluted at 2-8°C. DO NOT FREEZE! |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)