

Product datasheet for **RA023AP**

Feline IgG Feline Protein

Product data:

Product Type:	Native Proteins
Description:	Feline IgG feline protein, 1 mg
Species:	Feline
Protein Source:	Serum
Purity:	>95%
Buffer:	State: Liquid Buffer System: 0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4 containing 0.01% (w/v) Sodium Azide Presentation Label: AP
Preparation:	Liquid
Applications:	Suitable for antibody purification or to remove unwanted cross-reactivity.
Protein Description:	Feline IgG
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)