

Product datasheet for **RA023B**

Feline IgG Feline Protein

Product data:

Product Type:	Native Proteins
Description:	Feline IgG feline protein, 1 mg
Species:	Feline
Protein Source:	Serum
Concentration:	lot specific
Purity:	>95% prepared from normal serum delipidation, salt fractionation, ion exchange chromatography followed by extensive dialysis against the buffer stated above
Conjugation:	Biotin
Buffer:	State: Lyophilized Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) as stabilizer and 0.01% (w/v) Sodium Azide as preservative Label: Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Cat IgG molecule Presentation Label: Biotin
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent).
Preparation:	Lyophilized
Applications:	Immunoblotting. ELISA. Immunohistochemistry. Immunofluorescence.
Protein Description:	Feline IgG
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)