


Product datasheet for RA10107-150

P2X1 (P2RX1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens purinergic receptor P2X 1 (P2RX1)
Database Link:	Entrez Gene 25505 Rat
Synonyms:	P2X purinoceptor, Purinergic receptor, ATP receptor
Protein Families:	Druggable Genome, Ion Channels: ATP Receptors, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction

Product images:



[View online »](#)

