


## Product datasheet for RA10108-150

### P2X2 (P2RX2) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Monkey, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens purinergic receptor P2X 2 (P2RX2), transcript variant 6
Database Link:	<a href="#">Entrez Gene 114115 Rat</a> <a href="#">Entrez Gene 22953 Human</a>
Synonyms:	P2X purinoceptor 2, ATP receptor, Purinergic receptor
Protein Families:	Druggable Genome, Ion Channels: ATP Receptors, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction

#### Product images:



[View online »](#)

