


## Product datasheet for RA10109-150

### P2X3 (P2RX3) (C-term) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Monkey, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	purinergic receptor P2X 3
Database Link:	<a href="#">Entrez Gene 5024 Human</a>
Synonyms:	P2X purinoceptor 3, ATP receptor, Purinergic receptor

#### Product images:



DRG sections double-labeled for ASIC3 (red) and either P2X3 (green, A) or TRPV1 (green, B); double-labeled cells are yellow or orange



DRG sections double-labeled for ASIC3 (red) and either P2X3 (green, A) or TRPV1 (green, B); double-labeled cells are yellow or orange



[View online »](#)