

Product datasheet for RA103P

Avidin Chicken Protein

Product data:

Product Type:	Native Proteins
Description:	Avidin chicken protein, 10 mg
Species:	Chicken
Protein Source:	Egg white
Purity:	This product was prepared from chromatographically purified avidin.
Buffer:	State: Lyophilized purified Ig fraction without preservatives.
Bioactivity:	Specific: 14.4 U/mg protein (biotin titration method)
Reconstitution Method:	Restore with deionized water (or equivalent).
Preparation:	Lyophilized purified Ig fraction without preservatives.
Applications:	Suitable for conjugation or use as standard or control in many immunological assays. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Avidin. No reaction was observed against anti-Streptavidin.
Protein Description:	Unconjugated Avidin (Hen Egg White).
Storage:	Store vial at 2-8° C prior to restoration. For extended storage aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 2-8°C as an undiluted liquid. Dilute only prior to immediate use.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_990651
Locus ID:	396260
Summary:	avidin



[View online »](#)