


Product datasheet for RA12103-50

PGP9.5 (UCHL1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse, Porcine, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	ubiquitin C-terminal hydrolase L1
Database Link:	Entrez Gene 7345 Human
Synonyms:	UCH-L1, PGP 9.5, Ubiquitin thioesterase L1, Neuron cytoplasmic protein 9.5


Product images:




Endomorphin-2 Staining in rat dorsal horn.




[View online »](#)



END-2 staining of neurons from neonatal rat lumbar (L4/5) DRG. doi: 10.1111/j.1460-9568.2008.06238.x.



Endomorphin-2 Staining in rat dorsal horn.



END-2 staining of neurons from neonatal rat lumbar (L4/5) DRG. doi: 10.1111/j.1460-9568.2008.06238.x.