


## Product datasheet for RA14124-50

### Enkephalin (PENK) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens proenkephalin (PENK)
Database Link:	<a href="#">Entrez Gene 5179 Human</a>
Synonyms:	preproenkephalin; proenkephalin
Protein Families:	Secreted Protein


#### Product images:



Expression of ppENK in mouse striatal medium spiny neurons (MSNs). The Journal of Neuroscience, 8 September 2010, 30 (36): 11973-11982;doi: 10.1523/JNEUROSCI.2550-10.2010.



[View online »](#)



Expression of ppENK in mouse striatal medium spiny neurons (MSNs). *The Journal of Neuroscience*, 8 September 2010, 30 (36): 11973-11982;doi: 10.1523/JNEUROSCI.2550-10.2010.