

## Product datasheet for RA14139-100

### P2X3 (P2RX3) (C-term) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Primate, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	purinergic receptor P2X 3
Database Link:	<a href="#">Entrez Gene 5024 Human</a>
Synonyms:	P2X purinoceptor 3, ATP receptor, Purinergic receptor



[View online »](#)