


## Product datasheet for RA18009-100

### Synapsin I (SYN1) (pSer9) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens synapsin I (SYN1), transcript variant Ia
Database Link:	<a href="#">Entrez Gene 20964 MouseEntrez Gene 24949 RatEntrez Gene 6853 Human</a>
Synonyms:	Synapsin I, Brain protein 4.1, SYN1

#### Product images:



[View online »](#)

