


Product datasheet for RA18010-100

Synapsin I (SYN1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens synapsin I (SYN1), transcript variant la
Database Link:	Entrez Gene 20964 Mouse Entrez Gene 6853 Human
Synonyms:	Synapsin I, Brain protein 4.1, SYN1


Product images:



Synapsin and CREB staining of mouse retina.
Western blot analysis of mouse brain homogenates, untreated or phosphatase-treated, using Phospho-Synapsin (Ser9) (upper) or Synapsin Antibody (lower).



[View online »](#)



Synapsin and CREB staining of mouse retina. Western blot analysis of mouse brain homogenates, untreated or phosphatase-treated, using Phospho-Synapsin (Ser9) (upper) or Synapsin Antibody (lower).