


Product datasheet for RA18032-100

TCF4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IP, WB
Reactivity:	Chicken, Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	transcription factor 4
Database Link:	Entrez Gene 6925 Human
Synonyms:	Transcription factor 4, TCF-4, ITF2, ITF-2, SEF2, SEF-2


Product images:



Western blot analysis of extracts from various cell types using TCF4 (C48H11) Rabbit mAb.



[View online »](#)



Western blot analysis of extracts from various cell types using TCF4 (C48H11) Rabbit mAb.