


Product datasheet for RA18033-100

TCF4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IP, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	transcription factor 4
Database Link:	Entrez Gene 6925 Human
Synonyms:	Transcription factor 4, TCF-4, ITF2, ITF-2, SEF2, SEF-2


Product images:





TCF4 (C9B9) staining of HCT-116 cells (left) or SHSY-5Y cells (right) using TCF4 (C9B9) Rabbit mAb (green). Actin filaments have been labeled with Alexa Fluor® 555 phalloidin (red).




[View online »](#)



Western blot analysis of total cell lysates from HCT116, DLD1 and HepG2 cells using TCF4 (C9B9) Rabbit mAb.



TCF4 (C9B9) staining of HCT-116 cells (left) or SHSY-5Y cells (right) using TCF4 (C9B9) Rabbit mAb (green). Actin filaments have been labeled with Alexa Fluor® 555 phalloidin (red).



Western blot analysis of total cell lysates from HCT116, DLD1 and HepG2 cells using TCF4 (C9B9) Rabbit mAb.