


Product datasheet for RA19006-100

PARK7 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens Parkinsonism associated deglycase (PARK7), transcript variant 2
Database Link:	Entrez Gene 57320 Mouse Entrez Gene 11315 Human
Synonyms:	Oncogene DJ1, Parkinson disease protein 7
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Parkinson's disease

Product images:





Immunohistochemistry using anti-PARK7. Strong labeling of reactive astrocytes in Pick's disease (A) and adjacent to brain infarction (B).

PARK7 (DJ-1) Western Blot



[View online »](#)



Immunohistochemistry using anti-PARK7. Strong labeling of reactive astrocytes in Pick's disease (A) and adjacent to brain infarction (B).

PARK7 (DJ-1) Western Blot

