


Product datasheet for RA19012-100

IKAP (IKBKAP) (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase complex-associated protein (IKBKAP), transcript variant 1
Database Link:	Entrez Gene 8518 Human
Synonyms:	Elongator complex protein 1, IKAP, p150
Protein Families:	Druggable Genome


Product images:



Western blot of IKAP (1:500) in HEK293 cells stimulated with interleukin-1



[View online »](#)



Western blot of IKAP (1:500) in HEK293 cells stimulated with interleukin-1